


Arrow Peritoneal Lavage Kit



ARROW

Fast, simple, proven.



Suggested Procedure for Use

Note: This is an abbreviated pictorial Instruction for Use. Specific techniques and cautions apply. Please read product instructions for complete instructions.



Figure 1.

Insert 16 Ga. catheter/20 Ga. needle assembly through incision into peritoneum aiming toward the pelvic hollow. Advance tip of 16 Ga. catheter an additional 2-3 mm into peritoneal cavity.



Figure 2.

Remove 20 Ga. introducer needle.



Figure 3.

Insert desired end of spring-wire guide through 16 Ga. catheter into peritoneal cavity. Advance spring-wire guide into peritoneum to required depth.



Figure 4.

Remove introducer catheter and thread tip of 8 Fr. catheter over spring-wire guide (be certain that sufficient wire guide length remains exposed at hub end of catheter to maintain firm grip on wire guide). Grasping near skin, advance catheter to required depth in peritoneal cavity.

Hold catheter at required depth and remove spring-wire guide. Aspirate using 5cc syringe.



Figure 5.

Run solution into peritoneal cavity. Gently manipulate patient as required to insure proper intermixing of fluid.



Figure 6.

Lower bottlebag and I.V. tubing to floor level so that lavage fluid can drain freely by gravity due to siphoning effect.

Arrow

Peritoneal Lavage Kit

Arrow Peritoneal Lavage Kit allows diagnoses of intraperitoneal injury safely and simply by utilizing a percutaneous technique.

- Quickly performed and accurate diagnostic tool in the assessment of seriously injured patients with blunt abdominal trauma.¹
- Closed technique is tolerated well by the patient and is also extremely helpful in triage situations.
- Kit includes everything needed to easily and rapidly diagnose intraperitoneal bleeding with increased safety, reliability and convenience.
- J-tipped spring-wire guide for final catheter placement, making the older complication-prone trocar technique unnecessary.

Ordering Information

Product No.	Description	Per Case
AK-09000	Peritoneal Lavage Kit	10/case

One: Peritoneal Lavage Catheter: 8 Fr. x 9-1/4" (23.5 cm) Radiopaque FEP with Side Holes	Two: Syringe: 5 cc Luer-Slip
One: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)	One: 3 Swabstick Foil Pack 10% Povidone-Iodine
One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle	One: 5 ml Ampule HCL, 1% Lidocaine Solution
One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW	One: Drape: 18" x 26" with 3" fenestration
One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)	One: CSR Wrap
One: Injection Needle: 25 Ga. x 1" (2.54 cm)	Four: Gauze Pad(s): 4" x 4"
	Two: Gauze Pad(s): 2" x 2"
	One: Tubing Set: Fluid Administration
	One: Scalpel: #11

AK-09001	Peritoneal Lavage Kit	10/case
-----------------	------------------------------	----------------

One: Peritoneal Lavage Catheter: 8 Fr. x 9-1/4" (23.5 cm) Radiopaque FEP with Side Holes	One: Injection Needle: 25 Ga. x 1" (2.54 cm)
One: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)	Two: Syringe: 5 cc Luer-Slip
One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle	One: Drape: 18" x 26" with 3" fenestration
One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW	One: CSR Wrap
One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)	Four: Gauze Pad(s): 4" x 4"
	Two: Gauze Pad(s): 2" x 2"
	One: Tubing Set: Fluid Administration
	One: Scalpel: #11

For information regarding distribution in your area, contact:

ARROW
INTERNATIONAL

United States

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605
Phone: 610-378-0131
Call toll-free: 800-523-8446
Fax: 610-478-3199
Orders-only toll-free fax:
800-343-2935

Distribution Worldwide.

Offices in Africa, Canada, Czech Republic, France, Germany, Greece, Holland, Japan, Mexico, Slovakia and Spain.

www.arrowintl.com

For more information on Arrow's complete line of innovative products, call your Arrow representative or contact us by calling 800-523-8446 or 610-378-0131.

¹ Nagy KK, Krosner SM, Joseph KT, Roberts RR, Smith RF, Barrett J.

A method of determining peritoneal penetration in gunshot wounds to the abdomen. J Trauma. 1997;43:242-5

Caution: U.S. Federal Law limits this device to sale by or on order of a physician. Contents of unopened, undamaged package are sterile. Contains no medication. Disposable. Refer to package insert for current warnings, precautions, and instructions for use.

©2001 Arrow International, Inc. All rights reserved.

