

# Arrow Large Volume Paracentesis Tray

# ARROW

Alternative components

offer flow control choices.



Arrow  
Large Volume  
Paracentesis Tray

### Alternative components offer flow control choices.

- Rapid flow vacuum bottle tubing shortens procedure time while providing control
- Generous selection of components allows choice of drainage methods
- Polyurethane catheter with removable introducer needle and three-way stopcock assures patient safety and directed flow
- Comprehensive array of sterile components permits immediate use

### Ordering Information

Product No.	Description	Per Case
AK-00376	Large Volume Paracentesis Tray	6/case
One:	Vial: Vacutainer (EDTA)	Two: 5 ml Ampule HCl, 1% Lidocaine Solution
One:	Vial: Specimen	One: Dust Cap: Non-Vented
One:	Vial: Specimen	One: Drainage Bag
One:	Vial: Specimen	One: CSR Wrap
One:	Polyurethane Centesis Device: 8 Fr. Catheter over 18 Ga. x 7-1/2" (19 cm) RW Needle with 3-way Stopcock, Self-Sealing Valve and 5 cc Luer-Lock Syringe	One: Drape: 24" x 30" with 4" fenestration, with adhesive
One:	Introducer Needle: 20 Ga. x 3" (7.62 cm)	One: Dressing: Adhesive Strip 1 cm x 3 cm
Two:	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)	Five: Gauze Pad(s): 4" x 4"
One:	Injection Needle: 25 Ga. x 1" (2.54 cm)	One: Tubing Set: Drainage
One:	Syringe: 10 cc Luer-Lock	One: Scalpel: #11
One:	Syringe: 60 cc Luer-Lock	One: Tubing Set with Aspiration/Discharge Device
One:	3 Swabstick Foil Pack 10% Povidone-Iodine	Three: Label: Specimen Identification

For information regarding distribution in your area, contact:

**ARROW**  
INTERNATIONAL

**United States**  
Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605  
Phone: 610-378-0131  
Call toll-free: 800-523-8446  
Fax: 610-478-3199  
Orders-only toll-free fax:  
800-343-2935

*Distribution Worldwide.*  
*Offices in Africa, Canada, Czech Republic, France, Germany, Greece, Holland, Japan, Mexico, Slovakia and Spain.*

[www.arrowintl.com](http://www.arrowintl.com)

For more information on Arrow's complete line of innovative products, call your Arrow representative or contact us by calling 800-523-8446 or 610-378-0131.

Caution: U.S. Federal Law limits this device to sale by or on order of a physician. Contents of unopened, undamaged package are sterile. Contains no medication. Disposable. Refer to package insert for current warnings, precautions, and instructions for use.

©2001 Arrow International, Inc. All rights reserved.

PAR-C 09/01

