

Spring-Wire Guides

High pull strengths.

Arrow Spring-Wire Guides offer exceptional pull and tensile strength when compared to other brands. This strength helps assure more certain catheter placement with a minimum of risk.

In addition to high strength, Arrow was the first to introduce:

- Arrow DuoFlex® double distal with straight soft tips on both ends.
- “J” wire on one end with a straight soft tip on the other end.

Great selection.

Arrow offers a variety of Spring-Wire Guides, from .018" to .035" diameter, and lengths from 25 to 85cm.

A number of soft and “J” tips are offered.



Centimeter marked Spring-Wire Guides help determine insertion depth.



Two-piece Advancer and traditional straightener.



Arrow offers a broad range of sizes, diameters and tips.

High pull strength,
superior tensile strength
and a large variety to
meet your critical care
catheterization needs.

ARROW

Spring-Wire Guides

Ordering Information

| Product No. | Description | Per Case |
|-------------|--|----------|
| AW-01700 | .018" (.46mm) Dia. x 33-1/2" (85cm)—Straight Soft Tip on One End; "J" Tip on Other | 25 |
| AW-04018 | .018" (.46mm) Dia. x 9-3/4" (25cm)—Arrow DuoFlex® Straight Soft Tip on Both Ends | 25 |
| AW-04021-C | .021" (.53mm) Dia. x 13-13/16" (35cm)—Arrow DuoFlex® Straight Soft Tip on Both Ends | 25 |
| AW-04025 | .025" (.64mm) Dia. x 13-1/8" (33cm)—Arrow DuoFlex® Straight Soft Tip on Both Ends | 25 |
| AW-04025-J | .025" (.64mm) Dia. x 13-13/16" (35cm)—Straight Soft Tip on One End; "J" Tip on Other | 25 |
| AW-04032 | .032" (.81mm) Dia. x 12" (30.5cm)—Arrow DuoFlex® Straight Soft Tip on Both Ends | 25 |
| AW-04225 | .025" (.64mm) Dia. x 17-13/16" (45cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-04235 | .035" (.89mm) Dia. x 17-13/16" (45cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-04235-C | .035" (.89mm) Dia. x 17-13/16" (45cm)—Straight Soft Tip on One End; "J" Tip on Other | 25 |
| AW-04432 | .032" (.81mm) Dia. x 23-5/8" (60cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-04435 | .035" (.89mm) Dia. x 23-5/8" (60cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-14725 | .025" (.64mm) Dia. x 26-13/16" (68cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-14732 | .032" (.81mm) Dia. x 26-13/16" (68cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-14735 | .035" (.89mm) Dia. x 26-13/16" (68cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-16402 | .018" (.46mm) Dia. x 26-13/16" (68cm)—Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |
| AW-17752 | .035" (.89mm) Dia. x 39-1/2" (100cm)—PTFE-Coated Straight Soft Tip on One End; "J" Tip on Other—with Arrow Advancer™ | 25 |

PICC Wires

| Product No. | Description | |
|-------------|---|---|
| PW-18045 | .018" (.46mm) Dia. x 17-3/4" (45cm) Wire Proximal End: Solid Wire .018" Dia.; Distal End: Coil .018" Dia. x 9cm. | 5 |
| PW-18080 | .018" (.46mm) Dia. x 31-1/2" (80cm) Wire Proximal End: Solid Wire .018" Dia.; Distal End: Coil .018" Dia. x 7cm. | 5 |
| PW-18130 | .018" (.46mm) Dia. x 51-3/16" (130cm) Wire Proximal End: Solid Wire .018" Dia.; Distal End: Coil .018" Dia. x 7cm. | 5 |
| PW-18130-P | .018" (.46mm) Dia. x 51-3/16" (130cm) Wire Proximal End: Solid Wire .018" Dia.; Distal End: Coil .018" Dia. x 7cm; Coated with Platinum Tip | 5 |

Arterial Line Wires

| Product No. | Description | |
|-------------|---|----|
| RW-04018 | .018" (.46mm) Dia. x 5-3/8" (14cm) in Arterial Advancer Tube Assembly | 25 |
| RW-04025 | .025" (.64mm) Dia. x 5-3/8" (14cm) in Arterial Advancer Tube Assembly | 25 |

For information regarding distribution
in your area, contact:

INTERNATIONAL

United States
Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605
Phone: 610-378-0131
Call toll-free: 800-523-8446
Fax: 610-478-3199
Orders-only toll-free fax:
800-343-2935

Distribution Worldwide.
Offices in Africa, Canada, Czech Republic,
France, Germany, Greece, Holland, Japan,
Mexico, Slovakia and Spain.

www.arrowintl.com

Caution: U.S. Federal Law limits this device to sale by an order of a physician. Contents of unopened, undamaged package are sterile. Contains no medication. Disposable. Refer to package insert for current warnings, precautions, and instructions for use.

©2003 Arrow International, Inc. All rights reserved.
SWG-C 04/03

