

# Arrow TransRadial Artery Access Kit

*Now, access kits for physicians who prefer the easy-to-use Seldinger technique.*

**Quicker access because everything is included:**

- 21 Ga. introducer needle
- Arrow's new .018" floppy tip micro-puncture wire guide
- Choice of 5 Fr. or 6 Fr. introducer sheath sizes
- Sheath/dilator assembly

**Easier insertion for successful arterial access every time.**


- For easier sheath insertion, the dilator fits snugly over the micro-puncture wire guide, then gradually steps up to accommodate the OD of the sheaths.
- The sheaths, dilators and guidewires are precision-matched for consistent, successful arterial access.

**Smooth atraumatic insertion and removal without kinking or stretching.**

- ETFE thin-wall or new hydrophilically coated polyurethane sheath material promotes easy access.

**More options.**

- Catheter-Over-Needle Kiemeneij Technique Kits also available (see back page).



Everything you

need in one

complete kit.

Arrow TransRadial Artery Access Kit

Arrow  
TransRadial Artery  
Access Kit

## Seldinger Technique Kits with Micro-Puncture Wire Guide

Product No.	Contents	Packaged
AA-15611-S(new)	TransRadial Access Kit with 6 Fr. x 11cm Hydrophilically Coated Polyurethane Introducer Sheath and .018" Micro-Puncture Wire with 21 Ga. x 1-1/2" Introducer needle	5 kits/case
AA-15511	Arterial Access Kit with 5 Fr. x 11cm Introducer Sheath and .018" Micro-Puncture Wire with 21 Ga. x 1-1/2" Introducer needle	5 kits/case
AA-15611	TransRadial Access Kit with 6 Fr. x 11cm Introducer Sheath and .018" Micro-Puncture Wire with 21 Ga. x 1-1/2" Introducer needle	5 kits/case
AA-15624	TransRadial Access Kit with 6 Fr. x 24cm Introducer Sheath and .018" Micro-Puncture Wire with 21 Ga. x 1-1/2" Introducer needle	5 kits/case

## Micro-Puncture Wire Guides for Seldinger Technique

Product No.	Contents	Packaged
PW-18045	.018" x 45cm Straight Floppy Tip Wire Guide to introduce 11cm TransRadial Sheath/Dilator Assembly	5 /case
PW-18080	.018" x 80cm Straight Floppy Tip Wire Guide to introduce 24cm TransRadial Sheath/Dilator Assembly	5 /case

## Catheter-Over-Needle Technique Kits

Designed in conjunction with Dr. Ferdinand Kiemeneij – Interventional Cardiology-OLVG of Amsterdam, The Netherlands.

Product No.	Contents	Packaged
AA-00605(new)	TransRadial Access Kit with Arrow Radial Artery Catheterization System and 6 Fr. x 5cm Hydrophilically Coated Polyurethane Introducer Sheath with 150cm Exchange Wire Guide and Tegaderm® Dressing	5 kits/case
AA-00611	TransRadial Access Kit with Arrow Radial Artery Catheterization System and 6 Fr. x 11cm Introducer Sheath Assembly with 250cm Exchange Wire Guide	5 kits/case
AA-00624	TransRadial Access Kit with Arrow Radial Artery Catheterization System and 6 Fr. x 24cm Introducer Sheath Assembly with 250cm Exchange Wire Guide	5 kits/case

## Catheter-Over-Needle Technique Accessories

Product No.	Contents	Packaged
RA-04020	Arrow Radial Artery Catheterization Set with 20 Ga. x 1-3/4" FEP Radiopaque Dwelling Catheter over 22 Ga. Introducer Needle and Integral .018" Wire Guide, to access Radial Artery and inject Drug Cocktail	50/case
AV-04020	Arrow Arterial Catheterization Set with 20 Ga. x 1-3/4" FEP Radiopaque Dwelling Catheter over 22 Ga. Introducer Needle, to access Radial Artery and inject Drug Cocktail (20 Ga. Catheter will pass up to .025" Spring Wire Guide)	100/case

For information regarding distribution in your area, contact:

**ARROW**  
INTERNATIONAL

**United States**  
Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605  
Phone: 610-378-0131  
Call toll-free: 800-523-8446  
Fax: 610-478-3199  
Orders-only toll-free fax:  
800-343-2935

*Distribution Worldwide.*  
Offices in Africa, Canada, Czech Republic,  
France, Germany, Greece, Holland, Japan,  
Mexico, Slovakia and Spain.

[www.arrowintl.com](http://www.arrowintl.com)

For more information on Arrow's complete line of innovative products, call your Arrow representative or contact us by calling 800-523-8446 or 610-378-0131.

Caution: U.S. Federal Law limits this device to sale by or on order of a physician. Contents of unopened, undamaged package are sterile. Disposable. Refer to package insert for current warnings, precautions, and instructions for use.

\* Tegaderm is a registered trademark of 3M Company

©2001 Arrow International, Inc. All rights reserved.

TRA-C 04/01 11M

